

Lab Appendix A: Sharps and Biohazards - Safety Lancets and Finger Sticks

Sharps - Handling & Disposal

- Sharps are medical instruments that are used to puncture the skin (safety lancets*, lancets, needles, and blades)
- Dispose of sharps **immediately** in a hard red plastic sharps container!
- Do **NOT** dispose of sharps in a red biohazard bag

*Although safety lancets have retractable needles, they are still a safety hazard and must be disposed of in a sharps container.



Biohazard Bags - Handling & Disposal

- Bandages, used cotton and gauze, test strips, and gloves with body fluids on them are non-sharps biohazardous waste
- Dispose of non-sharps biohazardous waste in a red biohazard bag
- Wear gloves when handling biohazard bags



Safety Lancets

- Safety lancets have retractable needles or blades
- Are available with
 - Contact activation or
 - Triggered release
- Color coded by depth of penetration and gauge
 - Grey 1.8mm 23-Guage Low Flow 10-20
 - Orange 2.2mm 22-Guage Med Flow 20-40 ul
 - Pink 2.8mm 21-Guage Med-High Flow 40-60 ul



Surgilance® Safety Lancets

Where to Stick

Choose

- Non-dominant hand
- Inside the circles
- Black lines show lancet position

Best NOT to stick

- On the tip
- In the center
- Too far to the sides



Cleaning the Finger

- Use one firm swipe with an alcohol wipe (or isopropyl alcohol on a cotton ball)
 - Do not rub the alcohol wipe (or cotton ball) on the same area over and over again
- Let alcohol air dry – don't blow or fan
 - If alcohol is wet, the stick might sting
 - The sample should not be contaminated with alcohol
- Discard the alcohol wipe or cotton ball in a **biohazard bag**



Performing the Finger Stick

- Remove protective cap from the safety lancet



- Place red raised platform end of safety lancet on the test site



- Gently push the safety lancet down against test site to activate the safety lancet mechanism



- Once used, discard the safety lancet* into a **sharps container**

*Although safety lancets have retractable needles, they are still a safety hazard and must be disposed of in a sharps container.

